







UNIVERSAL BIPOLAR CABLE

2 in 1 Fits US and European bipolar instruments.
Reusable up to 100 autoclave cycles.

ORDERING INFORMATION

REF.	PRODUCT DESCRIPTION	PLUG		
P-30	Reusable cable for bipolar instruments. Plug for Erbe, Siemens Sutter BM50.		4 m	1 u.
Q-30	Reusable cable for bipolar instruments. Plug with 2 mm pins for Martin, Aesculap, Integra, Wolf, Berchtold, Sutter BM780.		4 m	1 u.
S-30	Reusable cable for bipolar instruments. Plug with twin 4 mm pin, Erbe, Eschmann, Alsa MCD series.		4 m	1 u.
V-30	Reusable cable for bipolar instruments. Block connector (28mm), 4 mm pins for Bowa, EMC, Erbe, Ethicon, Olympus, Soering, Conmed 5000.		4 m	1 u.

UNIVERSAL MONOPOLAR CABLE

2 in 1 4 / 8 mm Male plug & plug Female plug for 4-5 mm pin.
Reusable up to 100 autoclave cycles.



ORDERING INFORMATION

REF.	PRODUCT DESCRIPTION		
PMS8700	Monopolar Universal, 2 in 1 generator connector (8 mm / 4 mm) to Universal Female plug.	4 m	1 u.

This document is intended solely for the use of healthcare professionals according to his/her own professional clinical judgment. The information presented is only intended to demonstrate Prima's products, not to recommend a particular technique or procedure.

Products may not be available in all markets due to regulatory or commercial reasons.

No part of this document can be photocopied, reproduced or translated into any other languages without written consent.

Prima logo and products, are registered trademarks of Prima Medical Ltd. All other trademarks are trademarks of their respective owners or holders. The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Prima's trademark or other intellectual property rights concerning that name or logo. All rights reserved.