



BIPOLAR FORCEPS

- Lightweight and well balanced to prevent fatigue
- Highly polished tips minimise sticking
- Meets safety standard EN60601 requirements
- Non slip grip
- Fly lead with 3m cable
- Plug with 2 pins
- Single-use



ORDERING INFORMATION:


REF.	PRODUCT DESCRIPTION	
A FR-700-1105C	McPherson Bipolar Forceps, Straight, 0.5mm Tip, 10.8cm Length, 3m Cable, Single-Use	10
B FR-720-1210C	Adson Bipolar Forceps, 1.0mm Tip, 12.1cm Length, 3m Cable, Single-Use	10
C FR-730-1505C	Straight Bipolar Forceps, 0.5mm Tip, 15cm Length, 3m Cable, Single-Use	10
D FR-730-1510C	Straight Bipolar Forceps, 1.0mm Tip, 18cm Length, 3m Cable, Single-Use	10
D FR-730-1810C	Straight Bipolar Forceps, 1.0mm Tip, 20cm Length, 3m Cable, Single-Use	10
D FR-730-2010C	Straight Bipolar Forceps, 1.0mm Angled Tip, 20cm Length, 3m Cable, Single-Use	10
E FR-731-2010C	Straight Bipolar Forceps, 1.0mm Tip, 20cm Length, 3m Cable, Single-Use	10
F FR-750-1105C	Jeweler n.5 Bipolar Forceps, Straight, 0.5mm Tip, 11.4cm Length, 3m Cable, Single-Use	10
G FR-760-1610C	Jansen Bayonet Bipolar Forceps, 1.0mm Tip, 16.5cm Length, 3m Cable, Single-Use	10
G FR-760-2010C	Jansen Bayonet Bipolar Forceps, 1.0mm Tip, 20cm Length, 3m Cable, Single-Use	10
G FR-760-2210C	Jansen Bayonet Bipolar Forceps, 1.0mm Tip, 22cm Length, 3m Cable, Single-Use	10
G FR-760-2410C	Jansen Bayonet Bipolar Forceps, 1.0mm Tip, 24cm Length, 3m Cable, Single-Use	10

MONOPOLAR FORCEPS

- Well balanced to prevent hand fatigue
- Highly polished stainless-steel tips to minimise eschar build-up
- Meets safety standards
- Available in a pack of 10. Individually double - wrapped and sterile
- Pre-attached 3m cable with 8mm jack plug
- Fully insulated
- Single-use



ORDERING INFORMATION:

REF.	PRODUCT DESCRIPTION	
FR-380-1810C	Potts-Smith Monopolar forceps length 7" (17.8cm) 1.0mm tip serrated, 3m Cable, Single-Use	10
FR-380-1815C	Potts-Smith Monopolar forceps length 7" (17.8cm) 1.5mm tip serrated, 3m Cable, Single-Use	10
FR-380-2010C	Potts-Smith Monopolar forceps length 8" (20.3cm) 1.0mm tip serrated, 3m Cable, Single-Use	10
FR-380-2015C	Potts-Smith Monopolar forceps length 8" (20.3cm) 1.5mm tip serrated, 3m Cable, Single-Use	10
FR-380-2410C	Potts-Smith Monopolar forceps length 9.5" (24cm) 1.0mm tip serrated, 3m Cable, Single-Use	10
FR-380-2415C	Potts-Smith Monopolar forceps length 9.5" (24cm) 1.5mm tip serrated, 3m Cable, Single-Use	10
FR-390-2015C	Jansen Bayonet Monopolar forceps, length 8" (20.3cm) 1.5mm tip serrated, 3m Cable, Single-Use	10

This document is intended solely for the use of healthcare professionals according to his/her own professional clinical judgment. The information presented is only intended to demonstrate Prima's products, not to recommend a particular technique or procedure.

Products may not be available in all markets due to regulatory or commercial reasons.

No part of this document can be photocopied, reproduced or translated into any other languages without written consent.

Prima logo and products, are registered trademarks of Prima Medical Ltd. All other trademarks are trademarks of their respective owners or holders. The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Prima's trademark or other intellectual property rights concerning that name or logo. All rights reserved.